File Formats

Last update: 17 April 2018 Uview 15.4.0 and up

The following section describes the format of the proprietary file formats used in U-view. A simple program can be provided to list all the header info and LEEM overlay data in a table.

Still Image File

File Type: binary
File Extension: dat
Open with: U-view

Versions released:

Current FILEHEADER VERSION: 8

version 8: introduced with 3.3.0, 08/30/2010 CameraBitsPerPixel

version 7: introduced with 1.4.3, 03/21/06 recipe seperate: attachedRecipeSize

version 6: introduced with Uview 1.4.0, 11/29/05 image sequence

version 5: introduced with Uview 1.1.m, 2002 LEEMdata[256] added

Current IMAGEHEADER VERSION: 7

version 7: introduced 14.0.6, 2017 spareB2 changed to desired_rotation_fraction

Inclusion of MACROs in 'markup section' images are stored now without rotation but with infos desired rendering; desired rotation;

applied_processing

version 5: introduced with Uview 1.5.3, 03/19/2008 block of markeup for image tools

version 4: introduced with Uview 1.1.m, 2002: LEEMdata[256] added

Current LEEMDATAVERSION: 2

version 2: introduced with Uview 1.4.3, 03/21/06:

version 6: introduced with Uview 5.0, 6/22/2012

exposure time modified, average and seq. cycle added. *introduced with Uview 5.1, image header version 6:*

when leemdataversion>2 an extra block of LEEM data is added

File Structure overview:

- **1. FILE HEADER** (fixed size=104 Bytes)
- 2. **optional** for file version >=7: block of '**SEQUENCER RECIPE**' data (fixed size=128 Bytes) presence indicated by fileheader.attachedRecipeSize >0

3. IMAGEHEADER

- 4. **optional** for image version >=5: block of **IMAGE MARKUP** data (multiples of 128 Bytes) presence indicated by imageheader.attachedMarkupSize > 0
- 5. **optional** LEEM data if LEEMDataVersion>2 then size is LEEMDataVersion
- 6. **IMAGEDATA**: Width x Height x Pixel in 2 Bytes (if BitsPerPixels=16)

File structure (.dat & .dav) for file header version 7 and up:

Description Size file header FIXED: 104 Bytes only 1x per file + 128 Bytes if attachedRecipeSize >0 attached recipe else 0 Bytes (for.dav.always = 0) + image header FIXED: 288 Bytes if image header version >3 + markup block 0 Bytes if attachedMarkupSiz=0 X Bytes if attachedMarkupSiz>0 In that case X=128* ((attachedMarkupSize/128)+1) (integer division) + LEEM data 0 Bytes if *LEEMdataSize* else LEEMdataSize contains number of bytes in LEEM data block + image data ImageWidth*ImageHeight*2 Bytes

Notes for .dav files:

- 1. .day file consists of concatenated .dat files with the exception of the file header (only one)
- 2. attachedRecipeSize always=0
- 3. To get from file start to the first byte of the first image skip the following number of bytes:

104+ 288+ attachedMarkupSize+ LEEMdataSize

Offset from file start: attachedMarkupSize 126,127

LEEMdataSize 130,131

4. After reading the first image to get to the first byte if the following image skip the following number of bytes:

288+ attachedMarkupSize+ LEEMdataSize

Offset from file start: attachedMarkupSize 22,23

LEEMdataSize 26,27

byte offset from file start description

0 file start (start of file header)

.... start of image header #1

...

104+22 integer to calculate attachedMarkupSize0

• • •

104+26 integer containing LEEMdataSize0

. . .

begin markup block of size calculated from

 $attached Markup Size_0 \\$


...

104+288+ attachedMarkupSize begin LEEM data block of LEEMdataSize₀ length

. . .

104+288+ attachedMarkupSize₀ +**LEEMdataSize** begin image data of ImageHeight*ImageWith*2 length

. . .




C(++) program fragment to determine total size of (file) header for .dat with one image:

```
short s;
FILE *fd:
if((fd=fopen(currentfname, "r"))==NULL)return false;
fseek(fd,46,SEEK SET);
fread(&s,2,1,fd);
int attachedRecipeSize =(int)s;
fseek(fd,126 + attachedRecipeSize,SEEK_SET);
fread(&s,2,1,fd);
int attachedMarkupSize =(int)s;
attachedMarkupSize = 128*(( attachedMarkupSize /128)+1);
fseek(fd,130 + attachedRecipeSize ,SEEK SET);
fread(&s,2,1,fd);
int LEEMdataSize =(int)s;
int TotalSizeOfHeader=104+288+ attachedRecipeSize + attachedMarkupSize;
if(leemdatasize>2)TotalSizeOfHeader=TotalSizeOfHeader+ LEEMdataSize;
fclose(fd);
```

ConvertDATto ASCII.exe

This simple program opens .dat as well as .dav files and displays all header information. It also writes those information and optionally the data section as text into an output file.



Detailed description of data structure used in the still image file:

FILE HEADER

size in Bytes: 104		11
struct UKFileHeader{ char id[20]; short size; short version; short BitsPerPixel;	latest version bold 20 Bytes: contains "UKSOFT2001" + 0's 2 Bytes: sizeof(UKFileHeader) = 104 2 Bytes: FILEHEADER VERSION 2 Bytes, bits per pixel in file = 16	addressoffset 0 20 22 24
FILEHEADERVERSION>=8: short CameraBitsPerPixel; FILEHEADERVERSION<8: short spare;	2 Bytes, bits per pixel of camera	26
	s2 Bytes only !=0 for files taken by Uvie 1+1 Bytes only !=0 for files taken by Uvie 4 Bytes (because the following LONGLONG compiler at an 8-Byte boundary, this variable value of the file header descrip	w 8.5 and up 30 is aligned by the was not shown in
LONGLONG spare;	8 Bytes (was redundant 'starttime' in earl	y versions) 32
FILEHEADER VERSION>=2: short ImageWidth,ImageHeaper Short NrImages; FILEHEADER VERSION<2: short spare1,spare2,spare3;	eight; 2+2 Bytes 2 Bytes 6 Bytes	40,42 44
short attachedRecipeSize; attachedRecipeSize v =0: no a >0: attac	ttached block ched 128Byte block	<mark>46</mark>
BYTE spare[56];	ontains count of Bytes used for data within to 56 Bytes view 1.4.0,11/29/05 image sequence 'recipe of file) 2 Bytes 56 Bytes	48
short spare; BYTE spare[56];	2 Bytes 56 Bytes	104

SEQUENCER RECIPE block

Optional for file header version >=7

size in Bytes: 128

overview: sequencer commands arg1 arg2 arg3 <0 do nothing (<0 in first node: list is empty) 0 acquire image 1 wait msec 2 subtract image3 image= image2 >=Uview 6.1 3 calculate shift between image[arg1] & image[arg2] 4 set spin 1=spin up 0=spin down 5 normalize difference/sum 6 set LEEM supply absolute 7 set LEEM supply relative 8 LEEM preset 9 add and accumulate image= image2 image3 10 accumulate image[arg2] into internal buffer and display that buffer / nr cycles in arg 1 12 shift and accumulate image[arg2] into internal buffer and display that buffer in image[arg 1] Shift needs to be calculated prior to '12' >=Uview 6.1 filter-window 3x3... 13 apply median filter image >=Uview 6.1 14 divide image= factor*(image2 / image3) >=Uview 6.1

IMAGEHEADER
Image header version 5 & 6 & 7

size in Bytes: 288

struct UKImageHeader{	addressoff	
short size ;	2 Bytes: sizeof(UKImageHeader): 288	0
short version ;	2 Bytes: IMAGEHEADER VERSION	2
short ColorScaleLow, ColorScale		4,6
LONGLONG imagetime;	8 Bytes (see notes)	8
short MaskXShift, MaskYShift;	4 Bytes: shift of overlayed mask in x and y (pixels	
WORD RotateMask ;	2 Bytes: bit 7-15: 0-359 degrees image rotation,	20
	bits 0,1 UseMask	
	before Uview Version 1.6.5 not used (=0)	
short <mark>attachedMarkupSize</mark> ;	2 Bytes	22
attachedMarkupSize values:		
=0: no attached blo		
	KUP block is attached, value is number of	
	data within that block.	
·	his block on disk is a multiple of 128 bytes instead of just	
•	n previous versions. the actual length of the block calculate:	
	bytes: 128* ((attachedMarkupSize/128)+1)	
	arithmetic is used here i.e. 90/128=0)	
short spin ;	2 Bytes	24
short LEEMdataVersion ;	2 Bytes: = 1 or 2 or $>$ 2	26
	if >2 then extra block of LEEM data of the size eq	uals
	LEEMDataVersion is attached after the IMAGE	
	MARKUP block. This is implemented for	
	IMAGEHEADER VERSION>5.	
IMAGEHEADER VERSION>5		
IMAGEHEADER VERSION>5 unsigned char LEEMdata [239];	240 Bytes Overlay data	28
unsigned char LEEMdata[239];	240 Bytes Overlay data 1 Byte	28 267
	240 Bytes Overlay data 1 Byte 3 2 1 0	
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4	1 Byte	267
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4	1 Byte 3 2 1 0	267
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image	267
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values:	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2 LOG_RENDERING 3	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2 LOG_RENDERING 3 SQRT_RENDERING 4	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2 LOG_RENDERING 3 SQRT_RENDERING 4 ASinH_RENDERING 5	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2 LOG_RENDERING 3 SQRT_RENDERING 4 ASinH_RENDERING 5 GAUSS_RENDERING 6	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ; BYTE desired_rendering ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2 LOG_RENDERING 3 SQRT_RENDERING 4 ASinH_RENDERING 5 GAUSS_RENDERING 6 CLAHE_RENDERING 7	267 n s 268
unsigned char LEEMdata[239]; BYTE applied_processing; bit 5 4 bad pixel replacement background BYTE gray adjust zone; unsigned short backgroundvalue; BYTE desired_rendering;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2 LOG_RENDERING 3 SQRT_RENDERING 4 ASinH_RENDERING 5 GAUSS_RENDERING 6 CLAHE_RENDERING 7 fractional part of desired_rotation *100 version 7	267 n s 268
unsigned char LEEMdata [239]; BYTE applied_processing ; bit 5 4 bad pixel replacement background BYTE gray adjust zone ; unsigned short backgroundvalue ; BYTE desired_rendering ;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2 LOG_RENDERING 3 SQRT_RENDERING 4 ASinH_RENDERING 5 GAUSS_RENDERING 6 CLAHE_RENDERING 7 fractional part of desired_rotation *100 version 7 2 Bytes: additional argument for the selected	267 n s 268
unsigned char LEEMdata[239]; BYTE applied_processing; bit 5 4 bad pixel replacement background BYTE gray adjust zone; unsigned short backgroundvalue; BYTE desired_rendering; BYTE desired_rendering;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2 LOG_RENDERING 3 SQRT_RENDERING 4 ASinH_RENDERING 5 GAUSS_RENDERING 6 CLAHE_RENDERING 7 fractional part of desired_rotation *100 version 7 2 Bytes: additional argument for the selected rendering mode, not all modes use this argument	267 n s 268
unsigned char LEEMdata[239]; BYTE applied_processing; bit 5 4 bad pixel replacement background BYTE gray adjust zone; unsigned short backgroundvalue; BYTE desired_rendering;	1 Byte 3 2 1 0 subtracted reverseZ mirrorY mirrorX 90 degree rotation 1 Byte 8.6.4 0: invalid, 1=center half, 2= custom areas <0: entire image 2 Bytes 1 Byte Values: LINEAR_RENDERING 0 HISTogramEQUALisation_RENDERING 1 GAMMA_RENDERING 2 LOG_RENDERING 3 SQRT_RENDERING 4 ASinH_RENDERING 5 GAUSS_RENDERING 6 CLAHE_RENDERING 7 fractional part of desired_rotation *100 version 7 2 Bytes: additional argument for the selected	267 n s 268

```
short desired_rotation; 2 Bytes: angle the image was last displayed at short rotation_offset; 2 Bytes: for images<6 to avoid rotating them again short spare,spare; 4 Bytes

IMAGEHEADER VERSION=5
    unsigned char LEEMdata[256]; 256 Bytes Overlay data
    DWORD spare; 4 Bytes (to get to a total size divisible by 8)

notes:
```

The **'imagetime'** member of the image header structure holds the standard Windows *FILETIME* which is explained as the following:

"The *FILETIME* structure holds an unsigned 64-bit date and time value for a file. This value represents the number of 100-nanosecond units since the beginning of January 1, 1601." MS Visual C++ provides a number of functions to convert this time.

The **LEEMdata** array containing overlay data is structured as follows: source 1, argument 1 ... source n, argument n

source tags:

Highest bit of this byte is set when this item is <u>not</u> shown on the image: example: if 0x26 is tag for Start Voltage then 0xA6 is used as tag when Start voltage is recorded but not show on image.

```
0..99: 1. LEEM2000 module #
        2. followed by name
        3. followed by 1 ASCII digit identifying the unit
            unit codes:
            0=none,1=V,2=mA,3=A,4=C,5=K,6=mV,7=pA,8=nA,9=uA
        4. 0 to terminate the string
        5. data: 1 float (4 Bytes)
Oxff:
        skip
100:
        Mitutoyo micrometer readout: 2 floats (x, y coordinate)
101:
        (before 1.3.10) FOV (string max 16 char's + 0)
102:
        (before 1.3.10) varian controller #1 gauge #1 value (float)
103:
        (before 1.3.10) varian controller #1 gauge #2 value (float)
        camera exposure (float) in seconds - mislabeled previously as ms
104:
        LEEMdataVersion >1:
        2 bytes B1,B2 follow
           If B1>0 average is on , B2= number of averaged images 2 to 127
           if B1=0 average is off
           if B1<0 sliding average (<0 in this case -1: hex: 0xff, decimal:255)
105:
        title (string max 16 \text{ char's} + 0)
        varian controller #1 gauge #1 label, units value (string max 16
106:
                char's+0,+string max 4 char's+0, float) 27
107:
        varian controller #1 gauge #2 label, units value (string max 16
                char's+0,+string max 4 char's+0, float) 27
108:
        varian controller #2 gauge #1 label, units value (string max 16
                char's+0,+string max 4 char's+0, float) 27
109:
        varian controller #2 gauge #2 label, units value (string max 16
                char's+0,+string max 4 char's+0, float) 27
110:
        FOV, camera to FOV cal. factor (string max 16 char's+0+float)
111:
        phi, theta (float,float)
```

112:	spin	
113:	FOV rotation (from LEEM presets)	[15.2.0]
114:	Mirror state	[11.0.1]
115:	MCP screen voltage in kV	[11.0.1]
116:	MCP channelplate voltage in kV	[11.0.1]
120-1	30 additional gauges (#5,#6	

IMAGEHEADER (obsolete)

Image header version 4

Used in file header version 5 and 6 into file header version 7

size in Bytes: 288 struct UKImageHeader{ short size; 2 Bytes 2 Bytes short version; 4 Bytes filler LONGLONG imagetime; 8 Bytes BYTE spare[8]; 8 Bytes 2 Bytes short spin; short spareShort; 2 Bytes BYTE LEEMdata[256]; 256 Bytes 4 Bytes filler };

IMAGEHEADER (obsolete)

Image header version 3 and below Used in file header version <5

size in Bytes: 48

struct UKImageHeader{

short size; 2 Bytes short version; 2 Bytes

4 Bytes filler

LONGLONG imagetime; 8 Bytes long LEEMdata1_source; 4 Bytes float LEEMdata1_data; 4 Bytes short spin; 2 Bytes short spareShort; 2 Bytes float LEEMdata2_data; 4 Bytes BYTE spare[16]; 16 Bytes

};

IMAGE MARKUP block

Optional for file header version >=7

size in Bytes: 128, 256 etc. (multiples of 128)

size depends on actual space needed by the markers

The image header may be followed by 128 bytes or multiples of 128 bytes of an image markup block. This block contains info about lines and markers (letters & numbers) which the user has placed on the image.

The markup block consists of an array of 2-byte words:

[] description

- 0 block size
- 1 reserved
- 2 type of marker:
 - 1 horizontal cross section
 - 2 vertical cross section
 - 3 arbitrary cross section
 - 4 reserved
 - 6 markers
 - 7 inclusion or exclusion areas on image for histogram calculation (8.6.4)
 - 8 marker label (11.0.0) follows #6 immediately
 - 9 macro (from 14.0.0)

3... followed by:

in case of the cross sections:

First index x and y into image data array

Last index x and y into image data array

Cursor Data index x and y into cross section data

in case of the markers:

First index x and y into image data array

Last index x and y into image data array, y index ORed with type of markers

Color

Size and line weight

Followed in case of a marker containing text by the characters making up the text in case of inclusion or exclusion areas:

index of rectangle (0,1,2,3)

type: 1: inclusion, 2 exclusion

rectangle coordinates (left, top, right, bottom)

in case of marco:

followed by macro operation tag and arguments

note: the word 'markup' does not imply a similarity to html.

Data File containing multiple images

File Type: binary
File Extension: dat
Open with: U-view

File Contents:

Fileheader

Optional sequencer recipe block

Imageheader₁

Optional markup block₁
Optional LEEM data block₁

ImageData_{1:}

Width x Height x Pixel in 2 Bytes (if Bits PerPixels=16)

• • •

Imageheader_n

Optional markup data blockn Optional LEEM data blockn

ImageData_{1n}

Width x Height x Pixel in 2 Bytes (if Bits PerPixels=16)

Video File (.dav)

File Type: binary
File Extension: dav
Open with: U-view

File Contents: Concatenated still images (except that the file header appears only once)

Fileheader

Optional sequencer recipe block

Imageheader₁

Optional markup block₁ Optional LEEM data block₁

ImageData_{1:}

Width x Height x Pixel in 2 Bytes (if Bits PerPixels=16)

• • •

Imageheader_n

Optional markup data blockn Optional LEEM data blockn

ImageData_{1n}

Width x Height x Pixel in 2 Bytes (if Bits PerPixels=16)

Intensity data File

File Type: text File Extension: ivs

Open with: U-view or simple text editor : Microsoft Word-Pad

Notepad does not display the text correctly, Word does not save

correctly unless you follow the note below:

Note: don't re-save the file with a formatting editor like MS Word. If this is

done, the file can not be read back into U-view. You could use MS

number of data pairs to follow data pairs: time, intensity

Exponential format

Word but must save as plain text file.

File Contents:

UK SOFT

software FileVersion

IRectangle left top right bottom

StartChannel StartChannel

DataSection #Channels n

Time₁ intensity₁

....

Time_n intensity_n

last entry

example:

UK SOFT

software '

IRectangle 254 174 274 194

StartChannel 0 DataSection 4

5.050000e+003 1.251472e+006 5.220000e+003 1.252496e+006 5.270000e+003 1.253216e+006 5.380000e+003 1.254112e+006

last_entry